

taco™ Prep Bead Beater



- **Sample Capacity: 1 - 16**
- **Cycle duration: 13 sec/cycle & 40 sec/cycle**
- **Easy one-button operation**
- **Adjustable cycle number for effective homogenization**
- **Applicable to various biological samples**
- **Irregular collision of beads to achieve complete homogenization**
- **Fast & simple, no liquid nitrogen required**

taco™ Nucleic Acid Automatic Extraction System

A truly affordable extraction system



- **Cost Effective Magnetic Bead Purification**
- **Sample Capacity: 1 - 32 sample per run**
- **Easy one-button operation**
- **Universal built-in program**
- **Built in UV decontamination system**

Test Kits :

- **Taco DNA/RNA Extraction Kit 320 reaction**
- **Taco Plant DNA/RNA Extraction Kit, 320 reaction**
- **Taco Total DNA Extraction Kit, 320 reaction**